

# mpt64 Antibody



PACO52970

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mycobacterium tuberculosis

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:****Gene ID:**

mpt64

**Uniprot**

P9WIN9

**Synonyms:**

Immunogenic protein MPT64 (Antigen MPT64), mpt64

**Immunogen:**


Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: mpt64 antibody at 1.5 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 25 kDa. Observed band size: 25 kDa.