

ALAS1 Antibody, FITC conjugated



PACO52996

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, mitochondrial matrix, mitochondrion, nucleoplasm, 5-aminolevulinate synthase activity, identical protein binding, heme biosynthetic process, mitochondrion organization.

Gene ID:

ALAS1

Uniprot

P13196

Synonyms:

5-aminolevulinate synthase, nonspecific, mitochondrial (ALAS-H) (EC 2.3.1.37) (5-aminolevulinic acid, synthase 1) (Delta-ALA synthase 1) (Delta-aminolevulinate synthase 1), ALAS1, ALAS3 ALASH

Immunogen:

Recombinant Human 5-aminolevulinate synthase, nonspecific, mitochondrial protein (57-285AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A