

ETV1 Antibody, HRP conjugated



PACO53015

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional activator that binds to DNA sequences containing the consensus pentanucleotide 5'-CGGA[AT]-3'.

Gene ID:

ETV1

Uniprot

P50549

Synonyms:

ETS translocation variant 1 (Ets-related protein 81), ETV1, ER81

Immunogen:

Recombinant Human ETS translocation variant 1 protein (61-332AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A