

# ETV1 Antibody, Biotin conjugated



PACO53017

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcriptional activator that binds to DNA sequences containing the consensus pentanucleotide 5'-CGGA[AT]-3'.

**Gene ID:**

ETV1

**Uniprot**

P50549

**Synonyms:**

ETS translocation variant 1 (Ets-related protein 81), ETV1, ER81

**Immunogen:**

Recombinant Human ETS translocation variant 1 protein (61-332AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A