

TXNL4A Antibody, HRP conjugated



PACO53055

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Essential role in pre-mRNA splicing as component of the U5 snRNP and U4/U6-U5 tri-snRNP complexes that are involved in spliceosome assembly.

Gene ID:

TXNL4A

Uniprot

P83876

Synonyms:

Thioredoxin-like protein 4A (DIM1 protein homolog) (Spliceosomal U5 snRNP-specific 15 kDa protein) (Thioredoxin-like U5 snRNP protein U5-15kD), TXNL4A, DIM1 TXNL4

Immunogen:

Recombinant Human Thioredoxin-like protein 4A protein (1-142AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A