

# ZNF331 Antibody



PACO53066

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:200

**Protein Background:**

May be involved in transcriptional regulation. May play a role in spermatogenesis.

**Gene ID:**

ZNF331

**Uniprot**

Q9NQX6

**Synonyms:**

Zinc finger protein 331 (C2H2-like zinc finger protein rearranged in thyroid adenomas)  
(Zinc finger protein 361) (Zinc finger protein 463), ZNF331, RITA ZNF361 ZNF463

**Immunogen:**

Recombinant Human Zinc finger protein 331 protein (179-442AA).


**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images


---

N/A



N/A

Immunofluorescent analysis of HeLa cells using PACO53066 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human testis tissue using PACO53066 at dilution of 1:100.