

ZNF331 Antibody, FITC conjugated



PACO53068

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in transcriptional regulation. May play a role in spermatogenesis.

Gene ID:

ZNF331

Uniprot

Q9NQX6

Synonyms:

Zinc finger protein 331 (C2H2-like zinc finger protein rearranged in thyroid adenomas) (Zinc finger protein 361) (Zinc finger protein 463), ZNF331, RITA ZNF361 ZNF463

Immunogen:

Recombinant Human Zinc finger protein 331 protein (179-442AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A