

ECH1 Antibody, Biotin conjugated



PACO53101

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Isomerization of 3-trans,5-cis-dienoyl-CoA to 2-trans,4-trans-dienoyl-CoA.

Gene ID:

ECH1

Uniprot

Q13011

Synonyms:

Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial (EC 5.3.3), ECH1

Immunogen:

Recombinant Human Delta(3,5)-Delta(2,4)-dienoyl-CoA isomerase, mitochondrial protein (36-135AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A