

FAT3 Antibody



PACO53126

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May play a role in the interactions between neurites derived from specific subsets of neurons during development.

Gene ID:

FAT3

Uniprot

Q8TDW7

Synonyms:

Protocadherin Fat 3 (hFat3) (Cadherin family member 15) (FAT tumor suppressor homolog 3), FAT3, CDHF15 KIAA1989


Immunogen:

Recombinant Human Protocadherin Fat 3 protein (1209-1470AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO53126 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO53126 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).