

BIN3 Antibody, HRP conjugated



PACO53223

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in cytokinesis and septation where it has a role in the localization of F-actin.

Gene ID:

BIN3

Uniprot

Q9NQY0

Synonyms:

Bridging integrator 3, BIN3

Immunogen:

Recombinant Human Bridging integrator 3 protein (8-240AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A