

PTTG1IP Antibody



PACO53262

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May facilitate PTTG1 nuclear translocation.

Gene ID:

PTTG1IP

Uniprot

P53801

Synonyms:

Pituitary tumor-transforming gene 1 protein-interacting protein (Pituitary tumor-transforming gene protein-binding factor) (PBF) (PTTG-binding factor), PTTG1IP, C21orf1 C21orf3


Immunogen:

Recombinant Human Pituitary tumor-transforming gene 1 protein-interacting protein (118-180AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO53262 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO53262 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).