

GCLM Antibody, HRP conjugated



PACO53267

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, glutamate-cysteine ligase complex, enzyme regulator activity, glutamate-cysteine ligase catalytic subunit binding, glutamate metabolic process, glutathione biosynthetic process, positive regulation of glutamate-cysteine ligase activity, regulation of blood vessel size, response to drug, response to oxidative stress.

Gene ID:

GCLM

Uniprot

P48507

Synonyms:

Glutamate--cysteine ligase regulatory subunit (GCS light chain) (Gamma-ECS regulatory subunit) (Gamma-glutamylcysteine synthetase regulatory subunit) (Glutamate--cysteine ligase modifier subunit), GCLM, GLCLR

Immunogen:

Recombinant Human Glutamate--cysteine ligase regulatory subunit protein (1-274AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A