

RPS27A Antibody



PACO53334

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Gene ID:

RPS27A

Uniprot

P62979

Synonyms:

Ubiquitin-40S ribosomal protein S27a (Ubiquitin carboxyl extension protein 80)
[Cleaved into: Ubiquitin; 40S ribosomal protein S27a (Small ribosomal subunit protein eS31)], RPS27A, UBA80 UBCEP1


Immunogen:

Recombinant Human Ubiquitin-40S ribosomal protein S27a protein (16-56AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of MCF-7 cells using PACO53334 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).