

RPS27A Antibody, HRP conjugated



PACO53335

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	RPS27A
Source:	Uniprot
Rabbit	P62979
Isotype:	Synonyms:
IgG	Ubiquitin-40S ribosomal protein S27a (Ubiquitin carboxyl extension protein 80) [Cleaved into: Ubiquitin; 40S ribosomal protein S27a (Small ribosomal subunit protein eS31)], RPS27A, UBA80 UBCEP1
Applications:	Immunogen:
ELISA	Recombinant Human Ubiquitin-40S ribosomal protein S27a protein (16-56AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A