

# RPS27A Antibody, Biotin conjugated



PACO53337

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RPS27A

**Uniprot**

P62979

**Synonyms:**

Ubiquitin-40S ribosomal protein S27a (Ubiquitin carboxyl extension protein 80)  
[Cleaved into: Ubiquitin; 40S ribosomal protein S27a (Small ribosomal subunit protein eS31)], RPS27A, UBA80 UBCEP1

**Immunogen:**

Recombinant Human Ubiquitin-40S ribosomal protein S27a protein (16-56AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A