

PPIC Antibody, HRP conjugated



PACO53343

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

Gene ID:

PPIC

Uniprot

P45877

Synonyms:

Peptidyl-prolyl cis-trans isomerase C (PPIase C) (EC 5.2.1.8) (Cyclophilin C) (Rotamase C), PPIC, CYPC

Immunogen:

Recombinant Human Peptidyl-prolyl cis-trans isomerase C protein (26-188AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A