

RPL10A Antibody, Biotin conjugated



PACO53397

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, cytosol, cytosolic large ribosomal subunit, extracellular exosome, focal adhesion, membrane, nucleolus, nucleus, poly(A) RNA binding, RNA binding.

Gene ID:

RPL10A

Uniprot

P62906

Synonyms:

60S ribosomal protein L10a (CSA-19) (Large ribosomal subunit protein uL1) (Neural precursor cell expressed developmentally down-regulated protein 6) (NEDD-6), RPL10A, NEDD6

Immunogen:

Recombinant Human 60S ribosomal protein L10a protein (40-80AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A