

# TIE1 Antibody, Biotin conjugated



PACO53401

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transmembrane tyrosine-protein kinase that may modulate TEK/TIE2 activity and contribute to the regulation of angiogenesis.

**Gene ID:**

TIE1

**Uniprot**

P35590

**Synonyms:**

Tyrosine-protein kinase receptor Tie-1 (EC 2.7.10.1), TIE1, TIE

**Immunogen:**

Recombinant Human Tyrosine-protein kinase receptor Tie-1 protein (217-328AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A