

MSH2 Antibody, Biotin conjugated



PACO53421

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MSH2

Uniprot

P43246

Synonyms:

DNA mismatch repair protein Msh2 (hMSH2) (MutS protein homolog 2), MSH2

Immunogen:

Recombinant Human DNA mismatch repair protein Msh2 protein (2-254AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A