

BBOX1 Antibody, HRP conjugated



PACO53439

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the formation of L-carnitine from gamma-butyrobetaine.

Gene ID:

BBOX1

Uniprot

O75936

Synonyms:

Gamma-butyrobetaine dioxygenase (EC 1.14.11.1) (Gamma-butyrobetaine hydroxylase) (Gamma-BBH) (Gamma-butyrobetaine,2-oxoglutarate dioxygenase), BBOX1, BBH BBOX

Immunogen:

Recombinant Human Gamma-butyrobetaine dioxygenase protein (154-247AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A