

# BBOX1 Antibody, Biotin conjugated



PACO53441

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the formation of L-carnitine from gamma-butyrobetaine.

**Gene ID:**

BBOX1

**Uniprot**

O75936

**Synonyms:**

Gamma-butyrobetaine dioxygenase (EC 1.14.11.1) (Gamma-butyrobetaine hydroxylase) (Gamma-BBH) (Gamma-butyrobetaine,2-oxoglutarate dioxygenase), BBOX1, BBH BBOX

**Immunogen:**

Recombinant Human Gamma-butyrobetaine dioxygenase protein (154-247AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A