

PACO53446

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

Gene ID:

MIA

Uniprot

Q16674

Synonyms:

Melanoma-derived growth regulatory protein (Melanoma inhibitory activity protein), MIA


Immunogen:

Recombinant Human Melanoma-derived growth regulatory protein (87-124AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO53446 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate tissue using PACO53446 at dilution of 1:100.