

# MIA Antibody, Biotin conjugated



PACO53449

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

**Gene ID:**

MIA

**Uniprot**

Q16674

**Synonyms:**

Melanoma-derived growth regulatory protein (Melanoma inhibitory activity protein), MIA

**Immunogen:**

Recombinant Human Melanoma-derived growth regulatory protein (87-124AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A