

SCG3 Antibody, HRP conjugated



PACO53519

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasmic vesicle, extracellular region, secretory granule lumen, poly(A) RNA binding, platelet degranulation.

Gene ID:

SCG3

Uniprot

Q8WXD2

Synonyms:

Secretogranin-3 (Secretogranin III) (SgIII), SCG3

Immunogen:

Recombinant Human Secretogranin-3 protein (205-432AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A