

AK3 Antibody, Biotin conjugated



PACO53565

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in maintaining the homeostasis of cellular nucleotides by catalyzing the interconversion of nucleoside phosphates. Has GTP: AMP phosphotransferase and ITP: AMP phosphotransferase activities.

Gene ID:

AK3

Uniprot

Q9UIJ7

Synonyms:

GTP: AMP phosphotransferase AK3, mitochondrial (EC 2.7.4.10) (Adenylate kinase 3) (AK 3) (Adenylate kinase 3 alpha-like 1), AK3, AK3L1 AK6 AKL3L

Immunogen:

Recombinant Human GTP: AMP phosphotransferase AK3, mitochondrial protein (142-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A