

CCNE1 Antibody, FITC conjugated



PACO53620

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Essential for the control of the cell cycle at the G1/S (start) transition.

Gene ID:

CCNE1

Uniprot

P24864

Synonyms:

G1/S-specific cyclin-E1, CCNE1, CCNE

Immunogen:

Recombinant Human G1/S-specific cyclin-E1 protein (8-105AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A