

CDKN1C Antibody, FITC conjugated



PACO53636

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potent tight-binding inhibitor of several G1 cyclin/CDK complexes (cyclin E-CDK2, cyclin D2-CDK4, and cyclin A-CDK2) and, to lesser extent, of the mitotic cyclin B-CDC2.

Negative regulator of cell proliferation. May play a role in maintenance of the non-proliferative state throughout life.

Gene ID:

CDKN1C

Uniprot

P49918

Synonyms:

Cyclin-dependent kinase inhibitor 1C (Cyclin-dependent kinase inhibitor p57) (p57Kip2), CDKN1C, KIP2

Immunogen:

Recombinant Human Cyclin-dependent kinase inhibitor 1C protein (221-310AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A