

DCLK1 Antibody, Biotin conjugated



PACO53649

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable kinase that may be involved in a calcium-signaling pathway controlling neuronal migration in the developing brain. May also participate in functions of the mature nervous system.

Gene ID:

DCLK1

Uniprot

O15075

Synonyms:

Serine/threonine-protein kinase DCLK1 (EC 2.7.11.1) (Doublecortin domain-containing protein 3A) (Doublecortin-like and CAM kinase-like 1) (Doublecortin-like kinase 1), DCLK1, DCAMKL1 DCDC3A KIAA0369

Immunogen:

Recombinant Human Serine/threonine-protein kinase DCLK1 protein (143-381AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A