

ISL1 Antibody, Biotin conjugated



PACO53685

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds and regulates the promoters of the insulin, glucagon and somatostatin genes. Involved in the specification of motor neurons in cooperation with LHX3 and LDB1.

Gene ID:

ISL1

Uniprot

P61371

Synonyms:

Insulin gene enhancer protein ISL-1 (Islet-1), ISL1

Immunogen:

Recombinant Human Insulin gene enhancer protein ISL-1 protein (140-349AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A