

MEIS2 Antibody, FITC conjugated



PACO53696

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

MEIS2

Uniprot

O14770

Synonyms:

Homeobox protein Meis2 (Meis1-related protein 1), MEIS2, MRG1

Immunogen:

Recombinant Human Homeobox protein Meis2 protein (377-470AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A