

TRIM21 Antibody, HRP conjugated



PACO53755

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

TRIM21

Uniprot

P19474

Synonyms:

E3 ubiquitin-protein ligase TRIM21 (EC 2.3.2.27) (52 kDa Ro protein) (52 kDa ribonucleoprotein autoantigen Ro/SS-A) (RING finger protein 81) (RING-type E3 ubiquitin transferase TRIM21) (Ro(SS-A)) (Sjogren syndrome type A antigen) (SS-A) (Tripartite motif-containing protein 21), TRIM21, RNF81 RO52 SSA1

Immunogen:

Recombinant Human E3 ubiquitin-protein ligase TRIM21 protein (1-475AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A