

tal1 Antibody



PACO53778

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

tal1

Uniprot

O93507

Synonyms:

T-cell acute lymphocytic leukemia protein 1 homolog (TAL-1) (Stem cell protein) (zSCL), tal1, scl tal-1


Immunogen:

Recombinant Zebrafish T-cell acute lymphocytic leukemia protein 1 homolog protein (1-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Zebrafish tissue 20µg, 10µg. All lanes: tal1 antibody at 3µg/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 36, 23 kDa. Observed band size: 36 kDa.