

Vertebrate ancient opsin Antibody



PACO53826

Product Information

Size:

50ug

Reactivity:

Salmo salar

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

photoreceptor activity, phototransduction, visible light.

Gene ID:**Uniprot**

O13018

Synonyms:

Vertebrate ancient opsin


Immunogen:

Recombinant Salmo salar Vertebrate ancient opsin protein (1-75AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: Vertebrate ancient opsin antibody at 3 μ g/ml. Secondary, Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 25 kDa. Observed band size: 25 kDa.