

CYLD Antibody, Biotin conjugated



PACO53909

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CYLD

Uniprot

Q9NQC7

Synonyms:

Ubiquitin carboxyl-terminal hydrolase CYLD (EC 3.4.19.12) (Deubiquitinating enzyme CYLD) (Ubiquitin thioesterase CYLD) (Ubiquitin-specific-processing protease CYLD), CYLD, CYLD1 KIAA0849

Immunogen:

Recombinant Human Ubiquitin carboxyl-terminal hydrolase CYLD protein (269-435AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A