
Product Information

Size:

50ug

Reactivity:

Zebrafish, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

cytoplasm, type III intermediate filament.

Gene ID:

gfap

Uniprot

Q58EE9

Synonyms:

Glial fibrillary acid, c protein (GFAP), gfap


Immunogen:

Recombinant Zebrafish Glial fibrillary acid, c protein (1-444AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: MCF-7 whole cell lysate. All lanes: GFAP antibody at 3.4 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 50, 51 kDa. Observed band size: 50 kDa.