

MC4R Antibody, FITC conjugated



PACO54008

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor specific to the heptapeptide core common to adrenocorticotrophic hormone and alpha-, beta-, and gamma-MSH. Plays a central role in energy homeostasis and somatic growth. This receptor is mediated by G proteins that stimulate adenylate cyclase (cAMP).

Gene ID:

MC4R

Uniprot

P32245

Synonyms:

Melanocortin receptor 4 (MC4-R), MC4R

Immunogen:

Recombinant Human Melanocortin receptor 4 protein (1-43AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A