

NUMA1 Antibody, Biotin conjugated



PACO54033

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Highly abundant component of the nuclear matrix where it may serve a non-mitotic structural role, occupies the majority of the nuclear volume. Required for maintenance and establishment of the mitotic spindle poles during symmetric cell divisions, functioning as a tether linking bulk microtubules of the spindle to centrosomes. Also required for proper alignment of the mitotic spindle during asymmetric cell divisions.

Gene ID:

NUMA1

Uniprot

Q14980

Synonyms:

Nuclear mitotic apparatus protein 1 (Nuclear matrix protein-22) (NMP-22) (Nuclear mitotic apparatus protein) (NuMA protein) (SP-H antigen), NUMA1, NMP22 NUMA

Immunogen:

Recombinant Human Nuclear mitotic apparatus protein 1 protein (627-870AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A