

TM6SF2 Antibody, Biotin conjugated



PACO54065

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulator of liver fat metabolism influencing triglyceride secretion and hepatic lipid droplet content. May function as sterol isomerase.

Gene ID:

TM6SF2

Uniprot

Q9BZW4

Synonyms:

Transmembrane 6 superfamily member 2, TM6SF2

Immunogen:

Recombinant Human Transmembrane 6 superfamily member 2 protein (290-331AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A