

PACO54066

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

cytoplasm, intracellular membrane-bounded organelle, mitochondrion, nuclear membrane.

Gene ID:

TSPO

Uniprot

B1AH88

Synonyms:

Putative peripheral benzodiazepine receptor-related protein, TSPO, PBRS


Immunogen:

Recombinant Human Putative peripheral benzodiazepine receptor-related protein (12-85AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO54066 at dilution of 1:100.