

TSPO Antibody, FITC conjugated



PACO54068

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasm, intracellular membrane-bounded organelle, mitochondrion, nuclear membrane.

Gene ID:

TSPO

Uniprot

B1AH88

Synonyms:

Putative peripheral benzodiazepine receptor-related protein, TSPO, PBRS

Immunogen:

Recombinant Human Putative peripheral benzodiazepine receptor-related protein (12-85AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A