

TRBV7-9 Antibody, Biotin conjugated



PACO54101

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

plasma membrane, MHC protein binding, peptide antigen binding, immune response, regulation of immune response, T cell costimulation, T cell receptor signaling pathway.

Gene ID:

TCRB

Uniprot

P04435

Synonyms:

T-cell receptor beta chain V region CTL-L17, TCRB

Immunogen:

Recombinant Human T cell receptor β variable 7-9 protein (22-115AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A