

# AAGAB Antibody, HRP conjugated



PACO54135

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in endocytic recycling of growth factor receptors such as EGFR.

**Gene ID:**

AAGAB

**Uniprot**

Q6PD74

**Synonyms:**

Alpha- and gamma-adaptin-binding protein p34, AAGAB

**Immunogen:**

Recombinant Human Alpha- and gamma-adaptin-binding protein p34 protein (191-315AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A