

SEMA5B Antibody, FITC conjugated



PACO54192

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May act as positive axonal guidance cues.

Gene ID:

SEMA5B

Uniprot

Q9P283

Synonyms:

Semaphorin-5B, SEMA5B, KIAA1445 SEMAG

Immunogen:

Recombinant Human Semaphorin-5B protein (845-1013AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A