

BRD8 Antibody



PACO54206

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

BRD8

Uniprot

Q9H0E9

Synonyms:

Bromodomain-containing protein 8 (Skeletal muscle abundant protein) (Skeletal muscle abundant protein 2) (Thyroid hormone receptor coactivating protein of 120 kDa) (TrCP120) (p120), BRD8, SMAP SMAP2


Immunogen:

Recombinant Human Bromodomain-containing protein 8 protein (641-859AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry analysis of human lung cancer using PACO54206 at dilution of 1:100.