

AXIN2 Antibody, FITC conjugated



PACO54236

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	AXIN2
Source:	Uniprot
Rabbit	Q9BRF5
Isotype:	Synonyms:
IgG	Axin-2 (Axin-like protein) (Axil) (Axis inhibition protein 2) (Conductin), Axin2
Applications:	Immunogen:
ELISA	Recombinant Human AXIN2 protein (11-77AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A