

# MYRF Antibody, FITC conjugated



PACO54240

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Myelin regulatory factor: Constitutes a precursor of the transcription factor. Mediates the autocatalytic cleavage that releases the Myelin regulatory factor, N-terminal component that specifically activates transcription of central nervous system (CNS) myelin genes.

**Gene ID:**

MYRF

**Uniprot**

Q9Y2G1

**Synonyms:**

Myelin regulatory factor (EC 3.4) (Myelin gene regulatory factor) [Cleaved into: Myelin regulatory factor, N-terminal; Myelin regulatory factor, C-terminal], MYRF, C11orf9 KIAA0954 MRF

**Immunogen:**

Recombinant Human Myelin regulatory factor protein (843-1004AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A