

# KRT5 Antibody, HRP conjugated



PACO54431

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

cytoplasm, cytosol, extracellular exosome, intermediate filament, keratin filament, membrane, nucleus, scaffold protein binding, structural constituent of cytoskeleton, cornification.

**Gene ID:**

KRT5

**Uniprot**

P13647

**Synonyms:**

Keratin, type II cytoskeletal 5 (58 kDa cytokeratin) (Cytokeratin-5) (CK-5) (Keratin-5) (K5) (Type-II keratin Kb5), KRT5

**Immunogen:**

Recombinant Human Keratin, type II cytoskeletal 5 protein (316-491AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A