

MMP11 Antibody



PACO54434

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

May play an important role in the progression of epithelial malignancies.

Gene ID:

MMP11

Uniprot

P24347

Synonyms:

Stromelysin-3 (SL-3) (ST3) (EC 3.4.24) (Matrix metalloproteinase-11) (MMP-11), MMP11, STMY3


Immunogen:

Recombinant Human Stromelysin-3 protein (198-377AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of A549 cells using PACO54434 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).