

# MMP11 Antibody, Biotin conjugated



PACO54437

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May play an important role in the progression of epithelial malignancies.

**Gene ID:**

MMP11

**Uniprot**

P24347

**Synonyms:**

Stromelysin-3 (SL-3) (ST3) (EC 3.4.24) (Matrix metalloproteinase-11) (MMP-11), MMP11, STMY3

**Immunogen:**

Recombinant Human Stromelysin-3 protein (198-377AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A