

POLR2M Antibody



PACO54454

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:**Gene ID:**

POLR2M

Uniprot

Q6EEV4

Synonyms:

DNA-directed RNA polymerase II subunit GRINL1A, isoforms 4/5 (DNA-directed RNA polymerase II subunit M, isoforms 4/5), POLR2M, GRINL1A


Immunogen:

Recombinant Human DNA-directed RNA polymerase II subunit GRINL1A, isoforms 4/5 protein (2-91AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of A549 cells using PACO54454 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).